



21505

## Avidin ex. Egg White for molecular biology, 10U/mg Protein

Part B

CAS : 1405-69-2

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Activity	10 - 15 units/mg protein
Unit Definition	One unit is defined as the amount of avidin that will bind 1 µg D-biotin at pH 8.9.
Enzyme solution for activity determination	Prepared in 0.1 M phosphate buffer, pH 6.8.

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Mg	35021900 (GST 18%)
100 Mg	35021900 (GST 18%)
25 Mg	35021900 (GST 18%)

### Available Packages

10 Mg  
25 Mg  
100 Mg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com