



24766

## Arachis Oil pure

Part A

CAS : 68425-36-5

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Yellowish
Appearance (Form)	Viscous liquid
Density (g/ml) @ 20°C	~ 0.915 g/ml
Refractive Index (20°C)	1.468 - 1.470
Acid value	max. 0.5
Water	max. 0.3%

### Other Information

Notes	Suitable for the delivery of lipophilic compounds.
-------	--

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
250 ml	15089099 (GST 5%)
500 ml	15089099 (GST 5%)

### Available Packages

250 ml

500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com