



25903

## N-Methyl-2-Pyrrolidone (NMP) GC-HS, 99.9%

Part A

CAS : 872-50-4

Molecular Formula : C<sub>5</sub>H<sub>9</sub>NO

Molecular Weight : 99.13

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99.9 %
Residual solvents	
GC-HS Class 1	max. 1ppm
GC-HS Class 2	max. 10ppm
GC-HS Class 3	max. 50ppm
Density (g/ml) @ 20°C	1.029-1.031
Refractive Index (20°C)	1.470-1.471
Boiling Range	201-202°C
Non Volatile Matter	max. 0.0005%
Water (KF)	max. 0.05%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1000 ml	29142990 (GST 18%)
250 ml	29142990 (GST 18%)
500 ml	29142990 (GST 18%)
Type of Packing	
1000 ml	Glass Bottle
250 ml	Glass Bottle
500 ml	Glass Bottle

## Available Packages

250 ml

500 ml

1000 ml

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---