



26741

Ajowan Seed Oil extrapure

Part A

CAS : 8001-99-8

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Yellow to brown
Appearance (Form)	Liquid
Odour	Characteristic ajowan odour
Density (g/ml) @ 20°C	0.900 - 0.935 g/ml
Refractive Index (20°C)	1.488 - 1.505

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	33012590 (GST 18%)
500 ml	33012590 (GST 18%)

Available Packages

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com