



Sisco  
Research  
Laboratories

Product Datasheet  
and Specification

29986

## Pyruvate Kinase ex. Rabbit Muscle, 300U/mg powder

Part B

CAS : 9001-59-6

### Specifications

Appearance (Colour)	White
Appearance (Form)	Powder
Activity	min. 300 units /mg powder
Unit Definition	One unit will convert 1.0 micromole of phospho(enol)pyruvate to pyruvate per minute at pH 7.6 at 37°C.
Enzyme solution for activity determination	Prepared in potassium phosphate buffer, pH 7.6.

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
4000 Units	35079099 (GST 18%)

### Available Packages

2000 Units

4000 Units

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com