



Sisco  
Research  
Laboratories

Product Datasheet  
and Specification

35591

## Behenic Acid extrapure, 99%

Part B

CAS : 112-85-6

Molecular Formula : C<sub>22</sub>H<sub>44</sub>O<sub>2</sub>

Molecular Weight : 340.58

### Specifications

Appearance (Colour)	White to almost white
Appearance (Form)	Flakes or Powder or Beads or Chunks
Assay (GC)	min.99.0%
Melting Point	78-82°C

### General Information

Storage	8 to 25°C (Cool & Dry Area)
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Gms	38220090 (GST 12%)

### Available Packages

5 Gms  
25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com