



36678

## Orthophosphoric Acid extrapure, 85%

Part A

CAS : 7664-38-2

Molecular Formula : H<sub>3</sub>PO<sub>4</sub>

Molecular Weight : 98.00

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Density (g/ml) @ 20°C	1.70
Assay (T)	min. 85%
Chloride (Cl)	max. 0.002%
Sulphate (SO <sub>4</sub> )	max. 0.01%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.002%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1805 , IMCO Class No.:8 , Packing Group:III
HSN Code	
20 Ltr	28092010 (GST 18%)
2500 ml	28092010 (GST 18%)
500 ml	28092010 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

### Available Packages

500 ml

---

2500 ml

---

20 Ltr

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com