



51840

Capric Acid Methyl Ester (Methyl Decanoate) extrapure, 99%

Part A

CAS : 110-42-9

Molecular Formula : C₁₁H₂₂O₂

Molecular Weight : 186.29

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.871-0.873
Refractive Index (20°C)	1.424-1.425

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1000 ml	29159099 (GST 18%)
250 ml	29159099 (GST 18%)
Type of Packing	
250 ml	Glass Bottle
1000 ml	Glass Bottle

Available Packages

250 ml
1000 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com