



56874

## Tetramethylammonium Hydroxide 1M aq. solution ACS

Part A

CAS : 75-59-2

Molecular Formula : C<sub>4</sub>H<sub>13</sub>NO

Molecular Weight : 91.15

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Molarity	1.0M ±0.02M
Residue on ignition	max. 5ppm
Chloride (Cl)	max. 0.03%
Potassium (K)	max. 2ppm
Sodium (Na)	max. 2ppm

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1835 , IMCO Class No.:8 , Packing Group:III
HSN Code	
100 ml	38229090 (GST 12%)
250 ml	38229090 (GST 12%)
Type of Packing	
250 ml	Plastic Bottle

### Available Packages

100 ml

250 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com