



Sisco
Research
Laboratories

Product Datasheet
and Specification

59800

TCNT Type 1- Carbon Nanotubes Thin Walled, Length 5-30 μ m, OD 2-5nm, 90%

Part C

CAS : 308068-56-6

Specifications

Appearance (Colour)	Black
Appearance (Form)	Powder
Assay	min. 90%
OD	2 - 5 nm
Length	5 -30 μ m

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
100 Mg	28030090 (GST 18%)
25 Mg	28030090 (GST 18%)

Available Packages

25 Mg

100 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com