



60352

Alizarin Red S for tissue culture

Part B

CAS : 130-22-3

Molecular Formula : C₁₄H₇NaO₇S

Molecular Weight : 342.26

Specifications

Appearance (Colour)	Orange yellow
Appearance (Form)	Powder
Loss on drying	max. 8%
Absorbance (A) of 1% solution (pH 4.3) in a 1cm cell @424nm	.
pH range	pH 4.3 (brilliant yellow) to pH 6.3 (pink)
Suitability as indicator	Passes test
Suitability for Tissue Culture	Passes

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	32049000 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com