



60402

## L-Ascorbate Oxidase (Type 2) ex. Acremonium Sp., 350U/mg powder

Part B

CAS : 9029-44-1

### Specifications

Appearance (Colour)	Pale blue
Appearance (Form)	Powder
Activity	min. 350 units/mg powder
Unit Definition	One unit will oxidize 1.0 micromole of L-ascorbate to dehydroascorbate per min at pH 5.6 at 25 °C.
Enzyme solution for activity determination	Prepared in 0.1 M Phosphate/ EDTA buffer, pH 5.6

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1000 Units	35079099 (GST 18%)

### Available Packages

1000 Units

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com