



Sisco
Research
Laboratories

Product Datasheet
and Specification

64733

Behenic Acid pure, 85%

Part B

CAS : 112-85-6

Molecular Formula : C₂₂H₄₄O₂

Molecular Weight : 340.58

Specifications

Appearance (Colour)	White to almost white
Appearance (Form)	Flakes, Powder, Beads or Chunks
Purity by GC	min. 85%
Melting Point	74 - 78°C

General Information

Storage	8 to 25°C (Cool & Dry Area)
HSN Code	
25 Gms	38220090 (GST 12%)
250 Gms	38220090 (GST 12%)

Available Packages

25 Gms

250 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com