



66789

Protease ex. Bacillus Subtilis (Type 2 - Alkaline, granules), 100kAPU/g

Part B

CAS : 9036-06-0

Specifications

Appearance (Colour)	Off - white
Appearance (Form)	Granules
Activity	min. 100kAPU/g powder
Unit Definition	One unit will hydrolyze casein to produce color equivalent to 1.0 micromole (181 µg) of tyrosine per min at pH 7.5 at 37 °C (color by Folin-Ciocalteu reagent).
Enzyme solution for activity determination	Prepared in borax-NaOH buffer, pH 11.0.

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
250 Gms	35079099 (GST 18%)

Available Packages

250 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com