



81525

Protease ex. Bacillus Sp. (Type 5 - Neutral, granules), 8kNPU/g

Part B

CAS : 9036-06-0

Specifications

Appearance (Colour)	Off-white
Appearance (Form)	Granuals
Activity	min. 8kNPU/g powder
Unit Definition	One unit will hydrolyze casein to produce color equivalent to 1.0 micromole (181 µg) of tyrosine per min at pH 7.5 at 37 °C (color by Folin-Ciocalteu reagent).
Enzyme solution for activity determination	Prepared in borax-NaOH buffer, pH 11.0.

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	35079099 (GST 18%)
500 Gms	35079099 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

- 100 Gms
- 500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com